

Supplementary Table 1 Antibodies used

List of antibodies used for flow cytometry (FC), activation assays (AA), bead enrichment (BE) and immunohistochemistry (IHC)

<i>Target</i>	<i>Clone</i>	<i>Manufacturer</i>	<i>Application</i>
<i>CD3</i>	UCHT1	BD Biosciences	FC
<i>CD3</i>	OKT3	Invitrogen	AA
<i>CD28</i>	CD28.8	Invitrogen	AA
<i>CD45</i>	HI30	Tonbo Biosciences	FC
<i>CD45RA</i>	H100	BD Biosciences	FC
<i>CD45RO</i>	UCHL1	BD Biosciences	FC
<i>CD69</i>	FN50	BD Biosciences	FC
<i>CD103</i>	Ber-ACT 8	Miltenyi Biotec	FC
<i>TCR$\gamma\delta$ (pan)</i>	11F2	Miltenyi Biotec	FC
<i>TCRδ1</i>	REA173	Miltenyi Biotec	FC / BE
<i>TCRδ2</i>	B6	Biolegend	FC / BE
<i>TCRδ (pan)</i>	H-41	Santa Cruz Biotechnology	IHC
<i>IL17A</i>	eBio64DEC17	ThermoFisher	FC
<i>TNFα</i>	MAB11	BD Biosciences	FC
<i>CCR7</i>	150503	BD Biosciences	FC